

*100-31-181*

## CERTIFICATE OF ANALYSIS

Product:	Mouse Monoclonal Anti-Human CK-MM
Form:	Purified IgG from murine ascites
Catalog Number:	100-31-181
Clone Number:	MM1
Lot Number:	135MM1B
Purification:	Protein A affinity chromatography
Purity:	>98% by SDS-PAGE
Buffer:	Phosphate buffered saline, pH 7.2
Stabilizer:	None
Preservative:	0.05% sodium azide
Ig Class:	IgG <sub>2a</sub> ,k
Ig Concentration:	2.4 mg/mL
Long-Term Storage:	Frozen
Short-Term Storage:	Refrigerate
Shipping Temp:	Ambient

**\*\* FOR RESEARCH USE ONLY \*\***