

100-71-051

CERTIFICATE OF ANALYSIS

Product:	Mouse Monoclonal Anti-Human CD51
Form:	Purified IgG from murine ascites
Catalog Number:	100-71-051
Clone Number:	CD51.1
Lot Number:	32B35
Purification:	Protein A affinity chromatography
Purity:	>98% by SDS-PAGE
Buffer:	Phosphate buffered saline, pH 7.2
Stabilizer:	None
Preservative:	0.05% sodium azide
Ig Class:	IgG ₁ ,k
Ig Concentration:	2.6 mg/mL
Long-Term Storage:	Frozen
Short-Term Storage:	Refrigerate
Shipping Temp:	Ambient

**** FOR RESEARCH USE ONLY ****