



# PRODUCT DESCRIPTION

FOR IN VITRO RESEARCH USE ONLY  
NOT FOR USE IN HUMANS OR ANIMALS

ProteinTech Group, Inc  
2201 W. Campbell Park Dr.  
Chicago, IL 60612  
Toll Free: 1-888-4PTGLAB  
In Illinois: 1312-455-8498  
Fax: 1-312-455-8408  
Web: www.ptglab.com

## Antibody to STC2


**STCRP ; stanniocalcin 2**  
**Stanniocalcin-related protein**

**BC000658**

**About the antigen:**

### The antibody:

Format 1:	Affinity purified
Catalog No. :	10314-1-AP
Weight/Volume:	30 µg/150 µl
Tested Applications:	WB; IHC; IP




Human brain  
lysates were  
subjected to  
SDS PAGE  
followed by  
western blot with  
affinity purified Ab  
at dilution of  
1:1,500

### The production

Immunogen:	Hu, recombinant protein.
Tested Specificity:	Hu, Mo
Host:	Rb
Antibody species:	Polyclonal

### Storage Information:

The antibody is stored at 50% Glycerol in PBS buffer containing 0.02% sodium azide at -20 °C. The final concentration generally is 150 µg/ml for affinity purified antibody. Glycerol was added to a final concentration of 50% so that antibodies can be stored at -20 °C without going through freeze-thaw cycles during shipping.



Immunohistochemical staining of formalin-fixed paraffin-embedded breast cancer showing extracellular staining, with affinity purified Ab. Dilution :1:100