



# PRODUCT DESCRIPTION

FOR IN VITRO RESEARCH USE ONLY  
NOT FOR USE IN HUMANS OR ANIMALS

ProteinTech Group, Inc  
2201 W. Campbell Park Dr.  
Chicago, IL 60612  
Toll Free: 1-888-4PTGLAB  
In Illinois: 1312-455-8498  
Fax: 1-312-455-8408  
Web: www.ptglab.com

## DHCR24-BC004375 DHCR24

### BC004375

#### About the antigen:

#### The antibody:


Format 1: Affinity purified  
Catalog No. : 10471-1-AP  
Weight/Volume: 30 µg/150 µl  
Tested Applications: WB; IHC; IP

#### The production


Immunogen: Hu, recombinant protein.  
Tested Specificity: Hu, Mo  
Host: Rb  
Antibody species: Polyclonal

#### Storage Information:

The antibody is stored at 50% Glycerol in PBS buffer containing 0.02% sodium azide at -20 °C. The final concentration generally is 150 µg/ml for affinity purified antibody. Glycerol was added to a final concentration of 50% to so that antibodies can be stored at -20 °C without going through freeze-thaw cycles during shipping.



293 cell lysates were subjected to SDS PAGE followed by western blot with affinity purified Ab at dilution of 1:1,000



Immunohistochemical staining of formalin-fixed paraffin-embedded breast tumor showing cytoplasmic staining, with affinity purified Ab. Dilution :1:100