

Product Specification Sheet

Product: Non-Sterile Mouse Plasma in EDTA Potassium

Code: D408-05-0050

Lot # 11977

Size(s): 50 ml

Physical State: Liquid

Total Protein Concentration: 90 mg/ml (by Refractometry)

IgG Concentration: Normal

Anticoagulant Composition: This anticoagulant is composed of 1.5 mg/ml Potassium EDTA and 1.7 mg/ml Sodium Chloride, adjusted to a pH of 7.2.

Stabilizer: None

Preservative: None

Recommended Dilution(s): Not specified.

Storage Conditions: Store vial at -20° C prior to opening. Avoid cycles of freezing and thawing. Use standard sterile procedures when handling this product to maintain sterility. Expiration date is one (1) year from date of opening.

Sterilization: Not performed.

Immunoelectrophoresis: Normal

pH: Normal

Hemoglobin: Normal

Strain: Mixed breed and sex

Note: This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.