



Mouse Anti-Human CD95 (Fas/Apo-1)

Cat. No.	Form	Quantity
9735-01	Purified (UNLB)) Antibody	0.1 mg
9735-02	Fluorescein (FITC) Conjugate	100 tests
9735-08	Biotin (BIOT) Conjugate	100 tests
9735-09	R-phycoerythrin (R-PE) Conjugate	100 tests
9735-11	Allophycocyanin (APC) Conjugate	100 tests
9735-13	*Spectral Red™ (SPRD) Conjugate	100 tests
9735-14	Low Endotoxin, Azide-Free (LE/AF)	0.1 mg

DESCRIPTION

Clone	DX3
Ig Isotype	IgG _{2a}
Specificity	Human CD95/Fas/Apo-1

CD95, also known as Fas and Apo-1, is a 40-50 kDa type I transmembrane glycoprotein and a member of the tumor necrosis factor receptor superfamily. It is expressed by activated lymphocytes, monocytes, neutrophils, fibroblasts and cell lines. Fas ligand binding to CD95 induces apoptosis in activated mature lymphocytes thereby playing a role in maintaining peripheral tolerance. Crosslinking of CD95 by the monoclonal antibodies DX2 and DX3 delivers an apoptotic signal to Fas-sensitive cells, indicating that these monoclonal antibodies recognize a functional epitope of CD95.¹⁻⁵

RESEARCH APPLICATIONS

- Flow cytometry
- Immunohistochemistry (frozen sections)
- Immunoprecipitation
- *In vitro* induction of apoptosis

CHARACTERIZATION

To insure lot-to-lot consistency, each batch of product is tested by flow cytometry to conform with the characteristics of a standard reference.

WORKING DILUTIONS

Flow Cytometry:	Purified antibody	≤ 1 µg/10 ⁶ cells
	Fluorescein conjugate	10 µL/10 ⁶ cells
	Biotin conjugate	10 µL/10 ⁶ cells
	R-Phycoerythrin conjugate	10 µL/10 ⁶ cells
	Allophycocyanin conjugate	10 µL/10 ⁶ cells
	Spectral Red™ conjugate	10 µL/10 ⁶ cells


Other Applications: Since applications vary, you should determine the optimum working dilution of the product that is appropriate for your specific need.

For Research Use Only. Not for Diagnostic or Therapeutic Use.

IMMUNOFLUORESCENT STAINING

Product: Mouse Anti-Human CD95-FITC
Cat. No. 9735-02
Amount Used: 1 µg/10⁶ cells

The Jurkat human lymphoma cell line was stained with mouse anti-human CD95/Fas-FITC, following which large cells were gated and analyzed on a FACScan™ flow cytometer (BDIS, San Jose, CA).



HANDLING AND STORAGE

- The purified (UNLB) antibody is supplied as 1.0 mg of purified immunoglobulin in 1.0 mL of 100 mM borate buffered saline, pH 8.0. *No preservatives or amine-containing buffer salts added.* Store at 2-8°C.
- The fluorescein (FITC) conjugate is supplied as 100 tests in 1.0 mL PBS/NaN₃. Store at 2-8°C.
- The biotin (BIOT) conjugate is supplied as 100 tests in 1.0 mL of PBS/NaN₃. Store at 2-8°C.
- The R-phycoerythrin (R-PE) and allophycocyanin (APC) conjugates are supplied as 100 tests in 1.0 mL of PBS/NaN₃ and a stabilizing agent. Store at 2-8°C. **Do not freeze!**
- The Spectral Red™ (SPRD) conjugate is supplied as 100 tests in 1.0 mL of PBS/NaN₃ and a stabilizing agent. Store at 2-8°C. **Do not freeze!**
- The low endotoxin, azide-free (LE/AF) antibody is supplied as 0.5 mg purified immunoglobulin in 0.2 mL of PBS. **Aliquot and store at or below -20°C.**
- Conjugated forms should not be frozen and should be protected from prolonged exposure to light. Aliquot and freeze the low endotoxin, azide-free product at -20°C immediately upon receipt. Each reagent is stable for the period shown on the bottle label if stored as directed.

WARNING

Reagents contain sodium azide. Sodium azide is very toxic if ingested or inhaled. Avoid contact with skin, eyes, or clothing. Wear eye or face protection when handling. If skin or eye contact occurs, wash with copious amounts of water. If ingested or inhaled, contact a physician immediately. Sodium azide yields toxic hydrazoic acid under acidic conditions. Dilute azide-containing compounds in running water before discarding to avoid accumulation of potentially explosive deposits in lead or copper plumbing.

REFERENCES

1. Schlossman, S., L. Bloumsell, W. Gilks, J.M. Harlan, C. Kishimoto, J. Ritz, S. Shaw, R. Silverstein, T. Springer, T.F. Tedder, and R.F. Todd, eds. 1995. *Leukocyte Typing V: White Cell Differentiation Antigens*, Oxford University Press, Oxford.
2. Kishimoto, T., A.E.G. von dem Borne, S.M. Goyert, D.Y. Mason, M. Miyasaka, L. Moretta, K. Okumura, S. Shaw, T.A. Springer, K. Sugamura, and H. Zola, eds. 1998. *Leukocyte Typing VI: White Cell Differentiation Antigens*, Academic Press, New York.
3. Barclay, A.N., M.H. Brown, S.K.A. Law, A.J. McKnight, M.G. Tomlinson, and P.A. van der Merwe, eds. 1997. *The Leukocyte Antigens Facts Book, 2nd Edition*, CD95 Section, Academic Press, New York, p. 363.
4. Nagata, S., and P. Golstein. 1995. *Science* 267:3378.
5. van Parijs, L., and A.K. Abbas. 1996. *Curr. Opin. Immunol.* 8:355.

FACScan™ is a registered trademark of Becton Dickinson Immunocytometry Systems, Inc., San Jose, CA.

Spectral Red™ is a registered trademark of Southern Biotechnology Associates, Inc.

*Spectral Red™ is a PE/Cy™5 tandem conjugate. Cy™5 is a trademark of Amersham Biosciences Corp.. Cy™5 is for non-commercial research use only, not for therapeutic or in vivo applications. Other use needs license from Amersham Biosciences Corp., under U.S. Patent Nos. 4,981,977 and 5,268,486 and other patents pending. This material (or portions of this material) is subject to proprietary rights of Amersham Biosciences Corp. and Carnegie Mellon University and made and sold under license from Amersham Biosciences Corp.. This product is licensed for sale only for research. It is not licensed for any other use. There is no implied license hereunder for any commercial use. Commercial Use shall include; 1. sale, lease, license or other transfer of the material or any material derived or produced from it. 2. sale, lease, license or other grant of rights to use this material or any material derived or produced from it. 3. Use of this material to perform services for a fee for third parties. If you require a commercial license to use this material and do not have one, return this material, unopened to Southern Biotechnology, 160A Oxmoor Blvd, Birmingham, AL 35209, USA and any money paid for the material will be refunded.