

SUPPORT PRODUCTS

Control Immunoglobulins and Conjugates for Flow Cytometry

SANTA CRUZ BIOTECH

PRODUCT	IgG	FITC	PE	PE-Cy5	APC	PE-Cy7	APC-Cy7
normal mouse IgG ₁	sc-3877	sc-2855	sc-2866	sc-2877	sc-2888	sc-3770	sc-3781
normal mouse IgG _{2a}	sc-3878	sc-2856	sc-2867	sc-2878	sc-2889	sc-3771	sc-3782
normal mouse IgG _{2b}	sc-3879	sc-2857	sc-2868	sc-2879	sc-2890	sc-3772	sc-3783
normal mouse IgG ₃	sc-3880	sc-2858	sc-2869	sc-2880	sc-2891	sc-3773	sc-3784
normal mouse IgM	sc-3881	sc-2859	sc-2870	sc-2881	sc-2892	sc-3774	sc-3785
normal rat IgG ₁	sc-3882	sc-2827	sc-2871	sc-2882	sc-2893	sc-3775	sc-3786
normal rat IgG _{2a}	sc-3883	sc-2831	sc-2872	sc-2883	sc-2894	sc-3776	sc-3787
normal rat IgG _{2b}	sc-3884	sc-2835	sc-2873	sc-2884	sc-2895	sc-3777	sc-3788
normal rat IgM	sc-3885	sc-2839	sc-2874	sc-2885	sc-2896	sc-3778	sc-3789
normal Arm. hamster IgG	sc-3886	sc-2864	sc-2875	sc-2886	sc-2897	sc-3779	sc-3790
normal goat IgG	sc-3887	sc-3988	sc-3992	sc-3866	sc-3867	sc-3868	sc-3869
normal rabbit IgG	sc-3888	sc-3870	sc-3871	sc-3872	sc-3873	sc-3874	sc-3875

Affinity purified normal immunoglobulins and immunoglobulin conjugates for Flow Cytometry are provided for 100 tests for use as negative controls.