

## Phosphatase-Inhibitor-Mix I, powder

Mixture of 5 water-soluble inhibitors against acid and alkaline phosphatases, protein phosphatases 2A, 2B and 2C, phosphoprotein phosphatase, and protein-tyrosine phosphatase.

Contains imidazole, sodium fluoride, sodium molybdate, sodium-ortho-vanadate, and sodium tartrate.

The content of 1 vial dissolved in 1 ml water will give a 100-fold concentrate suitable for the treatment of 100 ml extract.

Storage Temperature: +15 °C to +30 °C

Cat.No.	Size
39050.01	1 vial
39050.02	5 vials
39050.03	10 vials