

JORDAN COMBS

COMBS

JSB-30 - COMBS

CATALOG#	DESCRIPTION	WELL WIDTH	SAMPLE VOLUME PER MM GEL
CSM-0808	0.8mm x 8 tooth	5.8mm	4.7 μ l
CSM-0813	0.8mm x 13 tooth	4.0mm	3.2 μ l
CSM-0818	0.8mm x 18 tooth	2.3mm	1.9 μ l
CSM-1508	1.5mm x 8 tooth	5.8mm	8.7 μ l
CSM-1513	1.5mm x 13 tooth	4.0mm	6.0 μ l
CSM-1518	1.5mm x 18 tooth	2.3mm	3.5 μ l

JSB-302 - COMBS

CATALOG#	DESCRIPTION	WELL WIDTH	SAMPLE VOLUME PER MM GEL
CSD-0816	0.8mm x 16 tooth	5.8mm	4.7 μ l
CSD-0827	0.8mm x 27 tooth	4.0mm	3.2 μ l
CSD-0836	0.8mm x 36 tooth	2.3mm	1.9 μ l
CSD-1516	1.5mm x 16 tooth	5.8mm	8.7 μ l
CSD-1522	1.5mm x 22 tooth	5.0mm	7.5 μ l
CSD-1527	1.5mm x 27 tooth	4.0mm	6.0 μ l
CSD-1533	1.5mm x 33 tooth	3.0mm	4.5 μ l
CSD-1536	1.5mm x 36 tooth	2.3mm	3.5 μ l

JSB-96 - COMBS

CATALOG#	DESCRIPTION	WELL WIDTH	SAMPLE VOLUME PER MM GEL
CB96-0825	0.8mm x 25 tooth, MCP**	5.8mm	4.7 μ l
CB96-0840	0.8mm x 40 tooth	4.0mm	3.2 μ l
CB96-0850	0.8mm x 50 tooth, MCP**	2.9mm	2.3 μ l
CB96-1525	1.5mm x 25 tooth, MCP**	5.8mm	8.7 μ l
CB96-1540	1.5mm x 40 tooth	4.0mm	6.0 μ l
CB96-1550	1.5mm x 50 tooth, MCP**	2.9mm	4.4 μ l
CB96-3025	3.0mm x 25 tooth, MCP**	5.8mm	17.4 μ l
CB96-3040	3.0mm x 40 tooth	4.0mm	12.0 μ l
CB96-15P	1.5mm, 2 markers, 1 sample	6.4, 204.0mm	9.6, 306.0 μ l
CB96-30P	3.0mm, 2 markers, 1 sample	6.4, 204.0mm	19.2, 612.0 μ l

JSB-75 - COMBS

CATALOG#	DESCRIPTION	WELL WIDTH	SAMPLE VOLUME PER MM GEL
CB75-0808	0.8mm x 8 tooth	9.9mm	7.9 μ l
CB75-0813	0.8mm x 13 tooth	5.8mm	4.7 μ l
CB75-0822	0.8mm x 22 tooth	4.0mm	16.0 μ l
CB75-0830	0.8mm x 30 tooth	2.5mm	2.0 μ l
CB75-1508	1.5mm x 8 tooth	9.9mm	14.9 μ l
CB75-1513	1.5mm x 13 tooth	5.8mm	8.7 μ l
CB75-1522	1.5mm x 22 tooth	4.0mm	6.0 μ l
CB75-1530	1.5mm x 30 tooth	2.5mm	3.7 μ l
CB75-3008	3.0mm x 8 tooth	9.9mm	29.72 μ l
CB75-3013	3.0mm x 13 tooth	5.8mm	17.5 μ l
CB75-3022	3.0mm x 22 tooth	4.0mm	12.0 μ l

JSB-120 - COMBS

CATALOG#	DESCRIPTION	WELL WIDTH	SAMPLE VOLUME PER MM GEL
CB120-0812	0.8mm x 12 tooth	10.0mm	8.1 μ l
CB120-0821	0.8mm x 21 tooth	5.8mm	4.7 μ l
CB120-0836	0.8mm x 36 tooth	4.0mm	3.2 μ l
CB120-0848	0.8mm x 48 tooth	2.5mm	2.0 μ l
CB120-1512	1.5mm x 12 tooth	10.0mm	15.2 μ l
CB120-1521	1.5mm x 21 tooth	5.8mm	8.7 μ l
CB120-1536	1.5mm x 36 tooth	4.0mm	6.0 μ l
CB120-1548	1.5mm x 48 tooth	2.5mm	3.7 μ l
CB120-3012	3.0mm x 12 tooth	10.0mm	30.5 μ l
CB120-3021	3.0mm x 21 tooth	5.8mm	17.5 μ l
CB120-3036	3.0mm x 36 tooth	4.0mm	12.0 μ l

MCP** - Multichannel Pipette Compatible

JVD-80 - COMBS

CATALOG#	DESCRIPTION	WELL WIDTH	VOLUME PER MM WELL DEPTH
CP-0406	0.4mm x 6 tooth	25.4mm	10.1 μ l
CP-0412	0.4mm x 12 tooth	9.9mm	4.0 μ l
CP-0418	0.4mm x 18 tooth	5.8mm	2.3 μ l
CP-0427	0.4mm x 27 tooth	4.0mm	1.6 μ l
CP-0806	0.8mm x 6 tooth	25.4mm	20.3 μ l
CP-0812	0.8mm x 12 tooth	9.9mm	7.9 μ l
CP-0818	0.8mm x 18 tooth	5.8mm	4.6 μ l
CP-0827	0.8mm x 27 tooth	4.0mm	3.2 μ l
CP-1506	1.5mm x 6 tooth	25.4mm	38.1 μ l
CP-1512	1.5mm x 12 tooth	9.9mm	14.9 μ l
CP-1518	1.5mm x 18 tooth	5.8mm	8.7 μ l
CP-1527	1.5mm x 27 tooth	4.0mm	6.0 μ l
CP-3006	3.0mm x 6 tooth	25.4mm	76.2 μ l
CP-3012	3.0mm x 12 tooth	9.9mm	29.7 μ l
CP-3018	3.0mm x 18 tooth	5.8mm	17.4 μ l
CP-3027	3.0mm x 27 tooth,	4.0mm	12.0 μ l
CP-04P	0.4mm, 2 markers, 1 sample	6.3, 152.4mm	2.5, 61.0 μ l
CP-08P	0.8mm, 2 markers, 1 sample	6.3, 152.4mm	5.0, 121.9 μ l
CP-15P	1.5mm, 2 markers, 1 sample	6.3, 152.4mm	9.5, 228.6 μ l
CP-30P	3.0mm, 2 markers, 1 sample	6.3, 152.4mm	18.9, 457.2 μ l
PF-26P	0.4mm x 26 tooth*	3.0mm	1.2 μ l
PF-36P	0.4mm x 32 tooth*	2.5mm	1.0 μ l
PF-46P	0.4mm x 46 tooth*	2.0mm	0.8 μ l

JQT-182 - COMBS

CATALOG#	DESCRIPTION	WELL WIDTH	SAMPLE VOLUME PER MM GEL
CQ-0418	0.4mm x 18 tooth	13.6mm	5.44 μ l
CQ-0432	0.4mm x 32 tooth	5.8mm	2.32 μ l
CQ-0454	0.4mm x 54 tooth	4.0mm	1.6 μ l
PF-64Q	0.4mm x 64 picket fence	1.9mm	0.8 μ l
PF-80Q	0.4mm x 80 picket fence	1.9mm	0.8 μ l



IBI ANALYTICAL COMBS



IBI PREPARATORY COMBS