

ACRYLIC PRODUCTS

ACRYLIC PRODUCTS
BETA PROTECTION

ORDERING INFORMATION

CAT.#	DESCRIPTION
KH-10	SMALL KIMWIPE™ HOLDER
KH-20	MEDIUM KIMWIPE™ HOLDER
KH-30	LARGE KIMWIPE™ HOLDER



KH-10, KH-20, & KH-30

PHYSICAL SPECIFICATIONS

KH-10

- **Outside Dimensions:** 6.25in. x 3.75in. x 6.0in.
- **Inside Dimensions:** 5.25in. x 3.0in. x 5.25in.
- **Kimwipe™ Towel:** 4.5in. x 8.375in.
- **Weight:** 1.5lbs

KH-20

- **Outside Dimensions:** 6.5in. x 4.5in. x 13.5in.
- **Inside Dimensions:** 5.5in. x 3.75in. x 12.75in.
- **Kimwipe™ Towel:** 11.75in. x 11.75in.
- **Weight:** 3.5lbs

KH-30

- **Outside Dimensions:** 6.5in. x 4.5in. x 16.5in.
- **Inside Dimensions:** 5.5in. x 3.75in. x 15.75in.
- **Kimwipe™ Towel:** 14.7in. x 16.6in.
- **Weight:** 4.5lbs

IBI's universal Kimwipe™ holders provide protection for your Kimwipe™ boxes on your bench top, while adding a bit of color to your lab. These versatile holders are supplied with an industrial strength double faced tape, mounting screws, and skid proof rubber feet; so mounting to the wall, sticking to cabinetry, or placing on your bench top, IBI Kimwipe™ holder is the product of choice for today's researcher!

Need Replacement Tape?
230230-00 - Replacement Tape - 5"
230234-00 - Replacement Tape - 9"



コスモ・バイオ株式会社