

# TUNAIR SHAKE FLASKS

MICROBIOLOGY  
PRODUCTS

## ORDERING INFORMATION

CATALOG#	DESCRIPTION
<b>SS-2011</b>	300ML NO-BAFFLE SHAKE FLASK-1EA.
<b>SS-4001</b>	300ML NO-BAFFLE SHAKE FLASK-4 PACK
<b>SS-2001C</b>	300ML NO-BAFFLE FLASK KIT (DRI-GAUZE LININGS)-1 KIT
<b>SS-2001S</b>	300ML NO-BAFFLE FLASK KIT (SILICONE LININGS)-1 KIT

*SS-2001C & SS-2001S flask kits come complete with twelve (12) polypropylene 300ml No-Baffle flasks, twelve (12) two piece caps, and five (15) filter linings.*

## WORKING VOLUMES

300ml Flasks . . . . . 100ml  
2.5L Flasks. . . . . 1L



## ORDERING INFORMATION

CATALOG#	DESCRIPTION
<b>SS-2012</b>	300ML HALF-BAFFLE SHAKE FLASK-1EA.
<b>SS-4002</b>	300ML HALF-BAFFLE SHAKE FLASK-4 PACK
<b>SS-2002C</b>	300ML HALF-BAFFLE FLASK KIT (DRI-GAUZE LININGS)-1 KIT
<b>SS-2002S</b>	300ML HALF-BAFFLE FLASK KIT (SILICONE LININGS)-1 KIT

*SS-2002C & SS-2002S flask kits come complete with twelve (12) polypropylene 300ml Half-Baffle flasks, twelve (12) two piece caps, and five (15) filter linings.*

## DIMENSIONS

Flask Size	Base Dia.	Neck Dia.	Height	Weight
300ml	3.25" [8.25cm]	1.75" [4.45cm]	6.00" [15.24cm]	.011b s [.004Kg]
2.5L	6.50" [16.51cm]	3.00" [7.62cm]	10.75" [27.31cm]	.021b s [.008Kg]

## ORDERING INFORMATION

CATALOG#	DESCRIPTION
<b>SS-2013</b>	300ML FULL-BAFFLE SHAKE FLASK-1EA.
<b>SS-4003</b>	300ML FULL-BAFFLE SHAKE FLASK-4 PACK
<b>SS-2003C</b>	300ML FULL-BAFFLE FLASK KIT (DRI-GAUZE LININGS)-1 KIT
<b>SS-2003S</b>	300ML FULL-BAFFLE FLASK KIT (SILICONE LININGS)-1 KIT

*SS-2003C & SS-2003S flask kits come complete with twelve (12) polypropylene 300ml Full-Baffle flasks, twelve (12) two piece caps, and five (15) filter linings.*

## CLEANING

All TUNAIR™ flasks and caps can be cleaned by soaking in water with a light detergent (e.g. Joy) solution to loosen dirt and contaminants, air dry. All TUNAIR™ flasks and caps can be autoclaved, if required.

## MIXING

No-Baffle (0 Baffles) . . . . . Normal Throws  
Half-Baffle (2 Baffles) . . . . . Vortex Motion  
Full-Baffle (6 Baffles) . . . . . Propeller Motion

## ORDERING INFORMATION

CATALOG#	DESCRIPTION
<b>SS-5011</b>	2.5L NO-BAFFLE SHAKE FLASK-1EA.
<b>SS-8001</b>	2.5L NO-BAFFLE SHAKE FLASK-4 PACK
<b>SS-6001C</b>	2.5L NO-BAFFLE FLASK KIT (DRI-GAUZE LININGS)-1 KIT
<b>SS-6001S</b>	2.5L NO-BAFFLE FLASK KIT (SILICONE LININGS)-1 KIT

*SS-6001C & SS-6001S flask kits come complete with four (4) polypropylene 2.5L No-Baffle flasks, four (4) two piece caps, and five (5) filter linings.*

## MATERIAL

All TUNAIR™ flasks and caps are constructed of chemical resistant polypropylene. All flasks and caps are fully autoclavable.



## ORDERING INFORMATION

CATALOG#	DESCRIPTION
<b>SS-5012</b>	2.5L HALF-BAFFLE SHAKE FLASK-1EA.
<b>SS-8002</b>	2.5L HALF-BAFFLE SHAKE FLASK-4 PACK
<b>SS-6002C</b>	2.5L HALF-BAFFLE FLASK KIT (DRI-GAUZE LININGS)-1 KIT
<b>SS-6002S</b>	2.5L HALF-BAFFLE FLASK KIT (SILICONE LININGS)-1 KIT

*SS-6002C & SS-6002S flask kits come complete with four (4) polypropylene 2.5L Half-Baffle flasks, four (4) two piece caps, and five (5) filter linings.*

## SHAKER SPEED

1" Throws . . . . . 150-200rpm  
2" Throws . . . . . 300-400rpm or possibly higher

## ORDERING INFORMATION

CATALOG#	DESCRIPTION
<b>SS-5013</b>	2.5L FULL-BAFFLE SHAKE FLASK-1EA.
<b>SS-8003</b>	2.5L FULL-BAFFLE SHAKE FLASK-4 PACK
<b>SS-6003C</b>	2.5L FULL-BAFFLE FLASK KIT (DRI-GAUZE LININGS)-1 KIT
<b>SS-6003S</b>	2.5L FULL-BAFFLE FLASK KIT (SILICONE LININGS)-1 KIT

*SS-6003C & SS-6003S flask kits come complete with four (4) 2.5L polypropylene Full-Baffle flasks, four (4) two piece caps, and five (5) filter linings.*