

TUNAIR SHAKE FLASKS

MICROBIOLOGY
PRODUCTS

ORDERING INFORMATION

| CATALOG# | DESCRIPTION |
|-----------------|---|
| SS-2011 | 300ML NO-BAFFLE SHAKE FLASK-1EA. |
| SS-4001 | 300ML NO-BAFFLE SHAKE FLASK-4 PACK |
| SS-2001C | 300ML NO-BAFFLE FLASK KIT (DRI-GAUZE LININGS)-1 KIT |
| SS-2001S | 300ML NO-BAFFLE FLASK KIT (SILICONE LININGS)-1 KIT |

SS-2001C & SS-2001S flask kits come complete with twelve (12) polypropylene 300ml No-Baffle flasks, twelve (12) two piece caps, and five (5) filter linings.

WORKING VOLUMES

300ml Flasks 100ml
2.5L Flasks. 1L



ORDERING INFORMATION

| CATALOG# | DESCRIPTION |
|-----------------|---|
| SS-2012 | 300ML HALF-BAFFLE SHAKE FLASK-1EA. |
| SS-4002 | 300ML HALF-BAFFLE SHAKE FLASK-4 PACK |
| SS-2002C | 300ML HALF-BAFFLE FLASK KIT (DRI-GAUZE LININGS)-1 KIT |
| SS-2002S | 300ML HALF-BAFFLE FLASK KIT (SILICONE LININGS)-1 KIT |

SS-2002C & SS-2002S flask kits come complete with twelve (12) polypropylene 300ml Half-Baffle flasks, twelve (12) two piece caps, and five (5) filter linings.

DIMENSIONS

| Flask Size | Base Dia. | Neck Dia. | Height | Weight |
|------------|--------------------|-------------------|---------------------|--------------------|
| 300ml | 3.25" [8.25cm] | 1.75" [4.45cm] | 6.00" [15.24cm] | .011bs [.004Kg] |
| 2.5L | 6.50" [16.51cm] | 3.00" [7.62cm] | 10.75" [27.31cm] | .021bs [.008Kg] |

ORDERING INFORMATION

| CATALOG# | DESCRIPTION |
|-----------------|---|
| SS-2013 | 300ML FULL-BAFFLE SHAKE FLASK-1EA. |
| SS-4003 | 300ML FULL-BAFFLE SHAKE FLASK-4 PACK |
| SS-2003C | 300ML FULL-BAFFLE FLASK KIT (DRI-GAUZE LININGS)-1 KIT |
| SS-2003S | 300ML FULL-BAFFLE FLASK KIT (SILICONE LININGS)-1 KIT |

SS-2003C & SS-2003S flask kits come complete with twelve (12) polypropylene 300ml Full-Baffle flasks, twelve (12) two piece caps, and five (5) filter linings.

CLEANING

All TUNAIR™ flasks and caps can be cleaned by soaking in water with a light detergent (e.g. Joy) solution to loosen dirt and contaminants, air dry. All TUNAIR™ flasks and caps can be autoclaved, if required.

MIXING

No-Baffle (0 Baffles) Normal Throws
Half-Baffle (2 Baffles) Vortex Motion
Full-Baffle (6 Baffles) Propeller Motion

ORDERING INFORMATION

| CATALOG# | DESCRIPTION |
|-----------------|--|
| SS-5011 | 2.5L NO-BAFFLE SHAKE FLASK-1EA. |
| SS-8001 | 2.5L NO-BAFFLE SHAKE FLASK-4 PACK |
| SS-6001C | 2.5L NO-BAFFLE FLASK KIT (DRI-GAUZE LININGS)-1 KIT |
| SS-6001S | 2.5L NO-BAFFLE FLASK KIT (SILICONE LININGS)-1 KIT |

SS-6001C & SS-6001S flask kits come complete with four (4) polypropylene 2.5L No-Baffle flasks, four (4) two piece caps, and five (5) filter linings.

MATERIAL

All TUNAIR™ flasks and caps are constructed of chemical resistant polypropylene. All flasks and caps are fully autoclavable.



ORDERING INFORMATION

| CATALOG# | DESCRIPTION |
|-----------------|--|
| SS-5012 | 2.5L HALF-BAFFLE SHAKE FLASK-1EA. |
| SS-8002 | 2.5L HALF-BAFFLE SHAKE FLASK-4 PACK |
| SS-6002C | 2.5L HALF-BAFFLE FLASK KIT (DRI-GAUZE LININGS)-1 KIT |
| SS-6002S | 2.5L HALF-BAFFLE FLASK KIT (SILICONE LININGS)-1 KIT |

SS-6002C & SS-6002S flask kits come complete with four (4) polypropylene 2.5L Half-Baffle flasks, four (4) two piece caps, and five (5) filter linings.

SHAKER SPEED

1" Throws 150-200rpm
2" Throws 300-400rpm or possibly higher

ORDERING INFORMATION

| CATALOG# | DESCRIPTION |
|-----------------|--|
| SS-5013 | 2.5L FULL-BAFFLE SHAKE FLASK-1EA. |
| SS-8003 | 2.5L FULL-BAFFLE SHAKE FLASK-4 PACK |
| SS-6003C | 2.5L FULL-BAFFLE FLASK KIT (DRI-GAUZE LININGS)-1 KIT |
| SS-6003S | 2.5L FULL-BAFFLE FLASK KIT (SILICONE LININGS)-1 KIT |

SS-6003C & SS-6003S flask kits come complete with four (4) 2.5L polypropylene Full-Baffle flasks, four (4) two piece caps, and five (5) filter linings.