

TUNAIR SHAKE FLASKS

MICROBIOLOGY
PRODUCTS

ORDERING INFORMATION

CATALOG#	DESCRIPTION
SS-2011	300ML NO-BAFFLE SHAKE FLASK-1EA.
SS-4001	300ML NO-BAFFLE SHAKE FLASK-4 PACK
SS-2001C	300ML NO-BAFFLE FLASK KIT (DRI-GAUZE LININGS)-1 KIT
SS-2001S	300ML NO-BAFFLE FLASK KIT (SILICONE LININGS)-1 KIT

SS-2001C & SS-2001S flask kits come complete with twelve (12) polypropylene 300ml No-Baffle flasks, twelve (12) two piece caps, and five (15) filter linings.

WORKING VOLUMES

300ml Flasks 100ml
2.5L Flasks. 1L



ORDERING INFORMATION

CATALOG#	DESCRIPTION
SS-2012	300ML HALF-BAFFLE SHAKE FLASK-1EA.
SS-4002	300ML HALF-BAFFLE SHAKE FLASK-4 PACK
SS-2002C	300ML HALF-BAFFLE FLASK KIT (DRI-GAUZE LININGS)-1 KIT
SS-2002S	300ML HALF-BAFFLE FLASK KIT (SILICONE LININGS)-1 KIT

SS-2002C & SS-2002S flask kits come complete with twelve (12) polypropylene 300ml Half-Baffle flasks, twelve (12) two piece caps, and five (15) filter linings.

DIMENSIONS

Flask Size	Base Dia.	Neck Dia.	Height	Weight
300ml	3.25" [8.25cm]	1.75" [4.45cm]	6.00" [15.24cm]	.011 b s [.004Kg]
2.5L	6.50" [16.51cm]	3.00" [7.62cm]	10.75" [27.31cm]	.021 b s [.008Kg]

ORDERING INFORMATION

CATALOG#	DESCRIPTION
SS-2013	300ML FULL-BAFFLE SHAKE FLASK-1EA.
SS-4003	300ML FULL-BAFFLE SHAKE FLASK-4 PACK
SS-2003C	300ML FULL-BAFFLE FLASK KIT (DRI-GAUZE LININGS)-1 KIT
SS-2003S	300ML FULL-BAFFLE FLASK KIT (SILICONE LININGS)-1 KIT

SS-2003C & SS-2003S flask kits come complete with twelve (12) polypropylene 300ml Full-Baffle flasks, twelve (12) two piece caps, and five (15) filter linings.

CLEANING

All TUNAIR™ flasks and caps can be cleaned by soaking in water with a light detergent (e.g. Joy) solution to loosen dirt and contaminants, air dry. All TUNAIR™ flasks and caps can be autoclaved, if required.

MIXING

No-Baffle (0 Baffles) Normal Throws
Half-Baffle (2 Baffles) Vortex Motion
Full-Baffle (6 Baffles) Propeller Motion

ORDERING INFORMATION

CATALOG#	DESCRIPTION
SS-5011	2.5L NO-BAFFLE SHAKE FLASK-1EA.
SS-8001	2.5L NO-BAFFLE SHAKE FLASK-4 PACK
SS-6001C	2.5L NO-BAFFLE FLASK KIT (DRI-GAUZE LININGS)-1 KIT
SS-6001S	2.5L NO-BAFFLE FLASK KIT (SILICONE LININGS)-1 KIT

SS-6001C & SS-6001S flask kits come complete with four (4) polypropylene 2.5L No-Baffle flasks, four (4) two piece caps, and five (5) filter linings.

MATERIAL

All TUNAIR™ flasks and caps are constructed of chemical resistant polypropylene. All flasks and caps are fully autoclavable.



ORDERING INFORMATION

CATALOG#	DESCRIPTION
SS-5012	2.5L HALF-BAFFLE SHAKE FLASK-1EA.
SS-8002	2.5L HALF-BAFFLE SHAKE FLASK-4 PACK
SS-6002C	2.5L HALF-BAFFLE FLASK KIT (DRI-GAUZE LININGS)-1 KIT
SS-6002S	2.5L HALF-BAFFLE FLASK KIT (SILICONE LININGS)-1 KIT

SS-6002C & SS-6002S flask kits come complete with four (4) polypropylene 2.5L Half-Baffle flasks, four (4) two piece caps, and five (5) filter linings.

SHAKER SPEED

1" Throws 150-200rpm
2" Throws 300-400rpm or possibly higher

ORDERING INFORMATION

CATALOG#	DESCRIPTION
SS-5013	2.5L FULL-BAFFLE SHAKE FLASK-1EA.
SS-8003	2.5L FULL-BAFFLE SHAKE FLASK-4 PACK
SS-6003C	2.5L FULL-BAFFLE FLASK KIT (DRI-GAUZE LININGS)-1 KIT
SS-6003S	2.5L FULL-BAFFLE FLASK KIT (SILICONE LININGS)-1 KIT

SS-6003C & SS-6003S flask kits come complete with four (4) 2.5L polypropylene Full-Baffle flasks, four (4) two piece caps, and five (5) filter linings.