



COSMO BIO CO., LTD.

Inspiration for Life Science


**Catalog Numbers**

SKB-IO-1

SKB-IO-2

## Ice On

TYPE 1 and TYPE 2



### Description

"Ice On" are aluminum-made racks for refrigeration of sample tubes. Samples can be kept cooler than directly inserting them into ice.

### Features and Advantages

Ice On Type 1 Cat#SKB-IO-1	Ice On Type 2 Cat#SKB-IO-2
<ul style="list-style-type: none"><li>• Tubes will not touch water and ice directly</li><li>• Made of aluminum and non-rust surface</li><li>• Compact design</li><li>• Can be stacked in a refrigerator or freezer.</li></ul>	
<ul style="list-style-type: none"><li>• Can cool 24 x 1.5 ml tubes</li><li>• Numbers (1-8) and alphabets (A to C) are inscribed on the surface. This makes the tubes easy to sort.</li></ul>	<ul style="list-style-type: none"><li>• Can cool 96 tubes of 0.2ml or 12x0.2ml 8-tube strips or 8x0.2ml 12-tube strips or 1x0.2ml 96 well plate and 6 tubes of 0.5ml and 4 tubes of 1.5ml.</li><li>• Easy to separate tubes of 1.5ml or 0.5ml with tubes of 0.2ml.</li></ul>

### Specifications

	Ice On Type 1 Cat#SKB-IO-1	Ice On Type 2 Cat#SKB-IO-2
<b>Material</b>	Aluminum dye cast	Aluminum dye cast
<b>Number of Samples</b>	24 tubes of 1.5ml (3 x 8)	96 tubes of 0.2ml tubes or 12 x 0.2ml 8 tube strips or 8 x 0.2ml 12 tube strips or 1 x 0.2ml 96 well plate 6 tubes of 0.5ml 4 tubes of 1.5ml
<b>Dimensions</b>	169 x 69 x 47mm	177 x 95 x 47mm
<b>Weight</b>	300g	470g



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

<http://www.cosmobio.com>

e-mail : [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

Phone : +81-3-5632-9617

FAX : +81-3-5632-9618