



PRODUCT CERTIFICATE

For: West Nile Virus Pre-M recombinant antigen

Catalog Number:	00402-V
Product Description:	West Nile Virus Pre-M (recombinant). N-terminal part of the Pre-M protein immunodominant region.
Purity:	~95% pure; Purity of proteins is evaluated by SDS-PAGE
Origin:	<i>E. coli</i>, recombinant
Fusions:	Six histidines.
Buffer/Composition:	20 mM Phosphate buffer pH 7.5.
Concentration:	1.0 mg/ml; bulk amount available upon request.
Application:	ELISA, WB, Lateral flow.
Storage:	Long Term (longer than three month): -80° C Short Term (up to three month): 4° C
Shipping:	Refrigerate upon arrival
Comments:	Immunoreactive with sera from WNV -infected individuals.

The products are for laboratory research use or further manufacturing only and should not be used for human therapeutic or diagnostic applications.

200 Dexter Avenue, Watertown, MA 02472 USA Phone: (617) 926-9167 Fax: (617) 926-9157
(800) VIROGEN (847-6436)

www.virogen.com