



PRODUCT CERTIFICATE

**For: HCMV (Human Cytomegalo Virus)
pp65 antibody**

Catalog Number:	023-A
Product Description:	HCMV pp65 IgG1 (monoclonal) Produced against a HCMV infected cells.
Specificity:	HCMV pp65 antigen.
Applications:	Western blot, ELISA, immunofluorecence.
Origin:	Mouse
Clone:	1-L-11
Purity:	Protein-A Sepharose purified.
Lot:	hcmv-pp65-7
Buffer/Composition:	1 X PBS pH 7.2; 0.01% sodium azide.
Concentration:	1.0 mg/ml; 100 µg/ vial, bulk amount available upon request.
Storage:	Long Term: -20°C Short Term: 4°C
Shipping:	Wet ice

The products are for laboratory research use or further manufacturing only and should not be used for human therapeutic or diagnostic applications.