



PRODUCT CERTIFICATE

For: Ad1 (Adenovirus) hexon E11 antibody

Catalog Number:	025-A
Product Description:	Ad1 hexon E11 IgG2a (monoclonal) Produced against purified Ad1 virus.
Specificity:	Ad1 hexon E11 antigen.
Applications:	Western blot (1:1,000), ELISA, Immunofluorescence.
Origin:	Mouse
Clone:	1-M-13
Purity:	Protein-A Sepharose purified.
Buffer/Composition:	1 X PBS pH 7.2; 0.01% sodium azide.
Concentration:	1mg/ml; 100 µg/ vial, bulk amount available upon request.
Storage:	Long Term: -20°C Short Term: 4°C
Shipping:	Wet ice

The products are for laboratory research use or further manufacturing only and should not be used for human therapeutic or diagnostic applications.