



PRODUCT CERTIFICATE

For: Ad7 (Adenovirus) hexon C11 antibody

Catalog Number:	026-A
Product Description:	Ad7 hexon C11 IgG2a (monoclonal) Produced against purified Ad7 virus.
Specificity:	Ad7 hexon C11 antigen.
Applications:	ELISA, Immunofluorecence.
Origin:	Mouse
Clone:	1-N-14
Purity:	Protein-A Sepharose purified.
Buffer/Composition:	1 X PBS pH 7.2; 0.01% sodium azide.
Concentration:	1mg/ml; 100 µg/ vial, bulk amount available upon request.
Storage:	Long Term: -20°C Short Term: 4°C
Shipping:	Wet ice

The products are for laboratory research use or further manufacturing only and should not be used for human therapeutic or diagnostic applications.