



PRODUCT CERTIFICATE

For: anti-FAS antibody

Catalog Number:	029-A
Product Description:	anti-FAS IgG1 (monoclonal) Produced against a recombinant protein.
Product Background:	Prevents apoptosis induced by actinomycin D and cycloheximide.
Specificity:	FAS antigen (human).
Origin:	Mouse
Clone:	6F-9
Purity:	Protein A purified.
Buffer/Composition:	1 X PBS pH 7.2; 0.01% sodium azide.
Concentration:	~ 1mg/ml; 100 µg/ vial, bulk amount available upon request.
Recommended Dilution:	For Western Blot: 1:1000 For ELISA: 1:1000
Storage:	Long Term: -20°C Short Term: 4°C
Shipping:	Wet ice
Application:	Immunoprecipitation (2-5 µg/10⁶ cells), Western Blotting (2-5 µg/ml), ELISA

The products are for laboratory research use or further manufacturing only and should not be used for human therapeutic or diagnostic applications.