



---

# PRODUCT CERTIFICATE

## For: anti-FAS Intracellular Fragment antibody

---

<b>Catalog Number:</b>	<b>039-A</b>
<b>Product Description:</b>	<b>anti-FAS IgG1 (monoclonal) Produced against a recombinant protein.</b>
<b>Specificity:</b>	<b>FAS antigen intracellular fragment (human).</b>
<b>Origin:</b>	<b>Mouse</b>
<b>Clone:</b>	<b>GH-9</b>
<b>Purity:</b>	<b>Protein A purified.</b>
<b>Buffer/Composition:</b>	<b>1 X PBS pH 7.2; 0.01% sodium azide.</b>
<b>Concentration:</b>	<b>~ 1mg/ml; 100 µg/ vial, bulk amount available upon request.</b>
<b>Recommended Dilution:</b>	<b>For Western Blot: 1:1000 For ELISA: 1:1000</b>
<b>Storage:</b>	<b>Long Term: -20°C Short Term: 4°C</b>
<b>Shipping:</b>	<b>Wet ice</b>
<b>Application:</b>	<b>Immunoprecipitation (2-5 µg/10<sup>6</sup> cells), Western Blotting (2-5 µg/ml), ELISA</b>

---

The products are for laboratory research use or further manufacturing only and should not be used for human therapeutic or diagnostic applications.