



PRODUCT CERTIFICATE

For: anti-FAS (CD95) monoclonal antibody

Catalog Number:	049-A
Product Description:	anti-FAS IgM (monoclonal)
Specificity:	Human FAS antigen (CD95). Induces apoptosis at 0.2 – 1.0 µg/ml.
Origin:	Mouse
Clone:	IPO-4
Applications:	Flowcytometry 0.1 - 1.0 µg/ml. Fluorescent microscopy 1.0 - 10.0 µg/ml. Immunohistology on frozen sections 1.0 – 10.0 µg/ml. Induction of Apoptosis 1.0 - 0.2 µg/ml.
Purity:	Protein G purified.
Buffer/Composition:	1 X PBS pH 7.2
Concentration:	1.0 mg/ml; 100 µg/ vial.
Storage:	Long Term: -20°C Short Term: 4°C
Shipping:	Wet ice

The products are for laboratory research use or further manufacturing only and should not be used for human therapeutic or diagnostic applications.