



---

## PRODUCT CERTIFICATE

### For: anti-CD45 RA monoclonal antibody

---

<b>Catalog Number:</b>	<b>067-A</b>
<b>Product Description:</b>	<b>anti-CD45 RA IgG2a (monoclonal)</b>
<b>Specificity:</b>	<b>CD45RA is an isoform, which consists of A, B, C determinants of CD45 (M<sub>r</sub> 220 kDa).</b>
<b>Origin:</b>	<b>Human</b>
<b>Clone:</b>	<b>158-4D3</b>
<b>Applications:</b>	<b>Flowcytometry. Fluorescent microscopy. Immunohistology on frozen sections. Western Blot.</b>
<b>Purity:</b>	<b>Protein G purified.</b>
<b>Buffer/Composition:</b>	<b>1 X PBS pH 7.2; 0.01%NaAzide.</b>
<b>Concentration:</b>	<b>1 mg/ml; 100 µg/ vial.</b>
<b>Storage:</b>	<b>Long Term: -20°C Short Term: 4°C</b>
<b>Shipping:</b>	<b>Wet ice</b>
<b>Refernces:</b>	<b>Schlossman S. F. et al. Leucocyte Typing V, Oxford University Press. 1995. p386-398</b>

---

The products are for laboratory research use or further manufacturing only and should not be used for human therapeutic or diagnostic applications.