



---

## PRODUCT CERTIFICATE

**For: anti-Nuclear Antigen 14-16 kDa (Proliferating Cells) antibody**

---

<b>Catalog Number:</b>	<b>079-A</b>
<b>Product Description:</b>	<b>anti-Nuclear Antigen (Proliferating Cells) IgM (monoclonal) Produced against a proliferating Harvey cell line.</b>
<b>Specificity:</b>	<b>Nuclear Antigen (human), 14-16 kDa.</b>
<b>Origin:</b>	<b>Mouse</b>
<b>Clone:</b>	<b>IPO-38</b>
<b>Purity:</b>	<b>Protein A purified.</b>
<b>Buffer/Composition:</b>	<b>1 X PBS pH 7.2</b>
<b>Concentration:</b>	<b>1mg/ml; 100 µg/ vial.</b>
<b>Recommended Dilution:</b>	<b>For ELISA: 1:1000</b>
<b>Storage:</b>	<b>Long Term: -20°C Short Term: 4°C</b>
<b>Shipping:</b>	<b>Wet ice</b>
<b>Application:</b>	<b>Immunoprecipitation, ELISA, Frozen Sections.</b>

---

The products are for laboratory research use or further manufacturing only and should not be used for human therapeutic or diagnostic applications.