

VI4

133A

anti-A1, A2, A3 blood antigen antibody

Catalog Number	133-A
Product Description	A1, A2, A3 IgM (monoclonal)
Clone	Z2B-1
Product Specificity	A1, A2, A3 human blood antigen.
Origin	Mouse hybridoma.
Purity	Protein A purified
Buffer Compositions	0.3 M NaCl; 1.25 mM EDTA; 0.1% Na Azide
Concentration	~1 mg/ml; 100 ug/vial, bulk amount available upon request
Application	Blood group typing in direct agglutination tests. - slide test - tube test - microplate test
Cross Reactivity	No cross-reactivity
Storage	Long Term: 4° C Short Term: 4° C Product is stable for 2 years at 4°C.
Expiration Date	Do not freeze
Shipping	Wet ice

These products are for laboratory research use or further manufacturing only and should not be used for human therapeutic or diagnostic applications.