



# PRODUCT CERTIFICATE

**For: anti-CD34 antibody**

---

<b>Catalog Number:</b>	<b>513-A</b>
<b>Product Description:</b>	<b>anti-CD34 IgG1 (monoclonal)</b>
<b>Specificity:</b>	<b>CD34 antigen (human). Hematopoietic precursor cells.</b>
<b>Origin:</b>	<b>Mouse</b>
<b>Clone:</b>	<b>ICO-115</b>
<b>Purity:</b>	<b>Protein A purified.</b>
<b>Buffer/Composition:</b>	<b>1 X PBS pH 7.4; 0.5% Saccharose; 0.2% gelatin; 0.1% NaN<sub>3</sub></b>
<b>Concentration:</b>	<b>1mg/ml; 100 µg/ vial.</b>
<b>Application:</b>	<b>Flow cytometry, IHC.</b>
<b>Storage:</b>	<b>4°C</b>
<b>Shipping:</b>	<b>Wet ice</b>

---

The products are for laboratory research use or further manufacturing only and should not be used for human therapeutic or diagnostic applications.