



PRODUCT CERTIFICATE

For: anti-CD38 antibody

| | |
|-----------------------------|---|
| Catalog Number: | 514-A |
| Product Description: | anti-CD38 IgG2a (monoclonal) |
| Specificity: | CD38 antigen (human). Hematopoietic precursors, thymocytes, activated T-cells, plasma cells. |
| Origin: | Mouse |
| Clone: | ICO-20 |
| Purity: | Protein A purified. |
| Buffer/Composition: | 1 X PBS pH 7.4; 0.5% Saccharose; 0.2% gelatin; 0.1% NaN₃ |
| Concentration: | 1mg/ml; 100 µg/ vial. |
| Application: | Flow cytometry. |
| Storage: | 4°C |
| Shipping: | Wet ice |

The products are for laboratory research use or further manufacturing only and should not be used for human therapeutic or diagnostic applications.