



VMRD, Inc.

Veterinary Medical Research and Development

Monoclonal Antibody

#BLV2

PO Box 502
Pullman, WA 99163 USA

Telephone Orders: (800) 221-8673
FAX Orders: (509) 331-5356
Technical Services: (509) 331-5815
<http://www.vmr.com>

Cell Line: BLV2
Specificity: Bovine Leukemia Virus
Isotype: IgG₁
Known Virus Reactivity: Anti-gp51; recognizes the D-D' epitope; sequential.
Known Applications: Radioimmunoprecipitation, Western blot, flow cytometry, antigen capture, immunoaffinity chromatography.

References:

No specific references at this time.

Description and Handling:

This monoclonal antibody is produced as mouse ascites fluid, clarified by centrifugation, and filtered through a 0.22 μ M filter. The concentration is 1 mg/ml in phosphate-buffered saline, preserved with sodium azide.

Storage:

When the vial is stored at 2-8°C, it should be stable for one year.

Warranty:

VMRD, Inc. warrants that the monoclonal antibodies are as described as to the quantity and contents stated on the label at the time of delivery to the customer. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE MADE BEYOND THE LABEL DESCRIPTION, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. Remedy is limited to replacement of the product or refund of the purchase price. VMRD, Inc. is not liable for property damage, personal injury, or economic loss caused by the product.

FOR RESEARCH PURPOSES ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC PROCEDURES.