

## Certificate of Analysis

Catalog No. COL2001/cell line: COLIS69A  
Specificity: Ecoli J5


Isotype: IgG1

Known applications: Flow cytometry, immunohistochemistry, cell function

Known species reactivity: Isotype control for IgG1

Description: The monoclonal antibody was produced in ascites fluid, clarified by filtration through a 0.2 micrometer filter. The antibody concentration is 1.0 mg/ml in phosphate buffered saline, preserved in 0.09% sodium azide.

Method used for quality control: Flow cytometry



Flow cytometric profile obtained with bovine leukocytes

**STORAGE:** When the vial is stored at 2-7°C, it should be stable for one year.

### REFERENCES:

Aydintug K, Inzana TJ, Letonja T, Davis WC, Corbeil LB. Cross-reactivity of monoclonal antibodies to *Escherichia coli* J5 with heterologous gram-negative bacteria and extracted lipopolysaccharides. *J. Infect.Dis.* 160:846-857, 1989.

### For research use only

Warranty: The WSU Monoclonal Antibody Center warrants this monoclonal antibody is as described and present in the amount stated on the label at the time of delivery to the customer. No other warranties, expressed or implied are made beyond the label description, including warranties of merchantability or fitness for a particular use. Remedy is limited to replacement of the product. The Center is not liable for property damage, personal injury, or economic loss caused by the product. The information listed in the data sheet is provided for reference only and not a substituted for the user's own respective material quality control.

