




---

**Monoclonal Antibody to human free kappa light chain**

**Code No.:** 7641  
**Quantity:** 200 $\mu$ l(200 $\mu$ g)  
**Price:** ¥40,000  
**Clone:** 5H11  
**Quality:** BALB/c mouse IgG2b $\kappa$  (Purified by protein A chromatography)  
**Solution:** 1mg/ml IgG, 0.1% NaN<sub>3</sub>, PBS, pH7.4  
**Character:**



**Storage:** The solution should be stored below -20°C.  
Repeated freezing and thawing should be avoided.

**Note:** Reasearch use only.