

Product Insert

TL (Tissue/cell lysate) Diluent Kit (96 Tests)

PN: P820301

COMPONENTS

TL Std Diluent. PN: 07TLS011. One bottle containing 2.5 mL of 1x standard diluent.

TL Assay Buffer. PN: 07TLA011. One bottle containing 8 mL of 1x assay buffer.

STORAGE: 2-8°C.

IMPORTANT: Sodium azide forms explosive compounds with heavy metals. These products contain <0.05% (w/w) azide which with repeated contact with lead and copper commonly found in plumbing drains may result in the buildup of shock sensitive compounds. Dispose in accordance with regulations from your institute.

APPLICATION: Optimal assay components for AimPlex immunoassays. To be used in conjunction with an AimPlex Immunoassay analyte kit and a basic kit. Refer to the AimPlex Multiplex Immunoassay User Manual and kit inserts for the assay procedure.

For Research Use Only. Not for use in diagnostic procedures.

Please contact Customer Service at 626-817-6366 with any questions or comments.

Thank you for your business!