

Product Insert

IGF I/IGF II Diluent Kit (96 Tests)

PN: P910102

COMPONENTS

IGF I/IGF II Stnd Diluent. PN: 07IS011.
One bottle containing 2.5 mL of 1x
standard diluent.

IGF I/IGF II Assay Buffer. PN: 07IA011.
One bottle containing 8 mL of 1x assay
buffer.

STORAGE: 2-8°C.

IMPORTANT: Sodium azide forms explosive compounds with heavy metals. These products contain <0.05% (w/w) azide which with repeated contact with lead and copper commonly found in plumbing drains may result in the buildup of shock sensitive compounds. Dispose in accordance with regulations from your institute.

APPLICATION: Optimal assay components for AimPlex Human and Non-Human Primate IGF I and IGF II 2-Plex Assay Panels. Applicable to all sample types listed in the assay kits. To be used in conjunction with an AimPlex Immunoassay analyte kit and a basic kit. Refer to the AimPlex Multiplex Immunoassay User Manual and kit inserts for the assay procedure.

For Research Use Only. Not for use in diagnostic procedures.

Please contact Customer Service at 626-817-6366 with any questions or comments.

Thank you for your business!